SERENITY NEXT

ISTRUZIONI PER L'USO
OPERATING INSTRUCTIONS
MODE D'EMPLOI
GEBRAUCHSANLEITUNG

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The manufacturer declines all and any responsibility for possible inexactitudes due to typing errors contained in this manual and reserves the right to modify his products without changing their main characteristics, if he considers such modifications necessary or useful.

To be preserved for future consultation.

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All polyurethane parts produced by Tecnodent are:

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The chair complies with EU Regulation no.
2017/745, concerning medical devices.

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INDEX

1.	Introduction	31
1.1.	Definitions.	31
1.2.	General warning.	31
1.3.	Documentation purpose.	31
1.4.	Document keeping.	31
1.5.	Responsibility limitations.	32
1.6.	Transport and storage.	32
1.7.	Packing disposal.	32
1.8.	Equipment disposal information.	32
2.	Equipment description	33
2.1.	Intended use.	33
2.2.	General features.	33
2.3.	Warranty.	34
2.4.	Identification and labelling.	34
2.5.	Technical specifications.	34
3.	Installation	37
3.1.	Unpacking.	37
3.2.	Supply status.	37
3.3.	Handling.	37
3.4.	Installation.	37
4.	Using the product	39
4.1.	Typical use.	39
4.2.	Controls.	39
4.3.	Movements.	40
5.	Maintenance	41
5.1.	General information.	41
5.2.	Cleaning.	41
5.3.	Main fuses replacement.	41
5.4.	Upholstery replacement.	42
6.	Troubleshooting	43
7.	Accessories	45
7.1.	Foot control with display.	45
7.2.	Protective paper towel holder.	46
7.3.	Upholstery with fire retardant and anti-static leatherette.	46
7.4.	Moulded memory foam upholstery	46
7.5.	Liquid containment basin, with valve.	46
7.6.	Electrically warmed seat upholstery.	47
7.7.	Back-up battery.	47
7.8.	Magnetic cushion.	47
7.9.	Integrated retractable leg rest.	48
7.10.	Patient side handles	48
7.11.	Side instruments holder.	48
7.12.	Transfusion / infusion armrest.	49
7.13.	Leg supports.	49
7.14.	Instrument tray holder support.	50
7.15.	Colposcope arm support.	51
7.16.	"Quick Release" system for backrest lowering.	51
7.17.	Additional electric socket on the base.	51
7.18.	Foot control cable wind up.	52
7.19.	Base with wheels.	52

1. INTRODUCTION

1.1. DEFINITIONS.

This manual adopts the following graphic and linguistic conventions.

NOTE.

Contains important information, which should be highlighted with respect to the text.

CAUTION.

Before conducting certain procedures, this message may be displayed. If not observed, it could cause damage to the equipment.

WARNING.

Before conducting certain procedures, this message may be displayed. If not observed, it could cause damage to the operator and the equipment.

Double function control(s).

Foot operated switches which control the multiple chair functions, in accordance to how they are activated. For example, continuous activation: manual chair movement; short activation (tap): automatic movement of a pre-programmed working position.

Synchronized movement.

Simultaneous and interdependent positioning of the chair parts, which may be activated in an automatic or manual setting (e.g., raising the chair and simultaneously leaning the backrest).

PCB.

Electronic card (printed circuit board).

CN.

Electric connector.

PWR.

Power supply.

1.2. GENERAL WARNING.

Upon user request, Tecnodent agrees to provide circuit diagrams, list of component parts, calibration instructions or other information available, which will assist the user's technical personnel to conduct repairs to those parts of the chair which are considered repairable by the manufacturer.

WARNING.

The manufacturer will accept responsibility for the safety, reliability and performance of the equipment on the following conditions:

- **The installation operations, modifications or repairs have been performed by authorised personnel;**
- **The electrical system of the installation site conforms to the legal regulations and standards;**
- **The equipment is utilized in accordance with the operating instructions.**

NOTE.

Please note that in accordance with Art. 14 of EEC Directive 85/374 "Responsibility for damages arising from defective products" implemented in Italy by Decreto del Presidente della Repubblica 24 maggio 1988, n. 224; "Right to compensation shall cease at the expiration of ten years from the day the manufacturer or the importer within the EU Nations has first marketed the product, which has caused the claimed damage".

1.3. DOCUMENTATION PURPOSE.

This manual shows the SERENITY NEXT gynecological chair and contains information regarding its correct use, functions, performances, maintenance, troubleshooting and corresponding solutions.

This manual is intended for the final user, the professional that uses the equipment to perform his job.

1.4. DOCUMENT KEEPING.

This manual is part of the equipment and must be conserved until its disposal.

It should be stored in a protected location, but near the device, for prompt usage if necessary; if this document is lost, the user should request a replacement copy.

1.5. RESPONSIBILITY LIMITATIONS.

The Manufacturer must be relieved of any responsibility in the following cases:

- Improper use or usage by untrained personnel for professional use;
- Use contrary to the applicable standards;
- Improper installation, if it is not specifically included in the supply contract;
- Power supply failure or defects;
- Lack of maintenance or servicing;
- Unauthorized modification or interventions;
- Use of non-original spare parts, or not model-specific;
- Failure to comply with the instructions (even partially);
- Exceptional events.

1.6. TRANSPORT AND STORAGE.

The equipment, if contained in its original packaging, may be exposed for a period of no more than 5 weeks to environmental conditions within the following ranges:

	Min	Max
Environment temperature:	-40° C	+70° C
Relative humidity (condensation included):	10%	95%
Atmospheric pressure:	500 hPa (mbar)	1060 hPa (mbar)

1.7. PACKING DISPOSAL.

All materials used for packaging are recyclable and respect the environment:

- Wooden pallet with fumigation treatment;
- Cardboard box;
- Polythene film or bubblewrap.

Please deliver the packing to an authorized landfill, for the recovery of these materials.

1.8. EQUIPMENT DISPOSAL INFORMATION.

According to the article 13 of D.lgs. no. 151, dated July 25th, 2005, "Implementation of the Directives 2002/95/EEC, 2002/96/EEC and 2003/108/EEC, concerning the reduction of hazardous substances in electric and electronic equipment, as well as waste disposal", and subsequent modifications such as D.lgs. no. 49, dated March 14, 2014, "Implementation of the Directive 2012/19/EEC concerning waste of electric and electronic equipment", Tecnodent states as follows.

- This symbol, placed in the identification label, indicates that the product at the end of its useful life must be disposed separately from other waste.
- The separate collection of end-of-life equipment is organized and managed by the manufacturer. The owner who wishes to dispose this equipment should contact the manufacturer or one of its representatives to allow separate collection of this end-of-life equipment.
- The correct separate collection for recycling, processing and environmentally friendly disposal helps avoiding potential adverse effects on human health, and promotes reuse and recycling of the equipment materials.
- Illegal disposal of the product by the owner will lead to penalties, as determined by the laws in force.
- If the product is placed out of use at the end of its working life, it is necessary to definitely put it out of service, by disconnecting the plug from the socket and cutting the power supply cable.

2. EQUIPMENT DESCRIPTION

2.1. INTENDED USE.

The product described by this manual is a patient chair for gynecological use. The equipment is designed for lifting and holding the patient in a proper working position during gynecological examinations and therapies.

It may be used by specialists in private clinics, hospital structures or specialized wards.

The apparatus requires the use of a compliant electrical installation, as specified in the section "installation". It is not advised for usage in close proximity with devices which monitor patient vital functions.

WARNING.

The equipment was developed and manufactured only for gynecological use. Any other use is prohibited. It is furthermore prohibited to perform any modification of the equipment or its parts, without explicit written authorization of the manufacturer.

NOTE.

The product has been designed to maintain the performance indicated in this manual during the service life, even in critical conditions specified by the standards. Any temporary deterioration (e.g., oscillations of the display in the case of electrical interference) will not compromise the device function.

2.2. GENERAL FEATURES.

2.2.1. Basic version.

A gynecological chair model SERENITY NEXT presents the following functions:

- Symmetric chair, for improved user versatility;
- Fixed base with levelling feet;
- Fully electro-mechanical movements, independent raising / lowering, backrest leaning, Trendelenburg and patient loading assistance (anti-Trendelenburg);
- Soft start/soft stop;
- Mechanical safety stop on the backrest;
- External paper roll holder, width 600 mm;
- Push-button hand-held panel;
- Moulded upholstery, free of seams;
- Leatherette with anti-microbial and anti-fungal properties;
- Low voltage controls (5 V);
- Removable liquid containment basin in stainless steel.

2.2.2. Accessories.

The chair can also be equipped with different accessories, depending on the working needs of the professional, as well as the comfort needs of the patient to be treated. In detail, the following optionals are available:

- Operator foot rest;
- Side instruments holders;
- Base with wheels (for short distance movements e.g. during cleaning operations);
- Foot control with display;
- Memorized work positions (n°3+1 Program x n° 3 profiles), available only with foot control;
- Automatic return to the patient exit position, available only with foot control;
- "SLOW MODE" setting, available only with foot control;
- Fire retardant and antistatic upholstery;
- Magnetic cushion;
- Patient side handles (movable);
- Integrated extractable leg rest, for examination bed use;
- Removable leg supports, with under thigh (Göpel) holder;
- Removable leg supports, with under feet hoolder;
- Motorised leg supports, with under feet hoolder;
- "Quick Release" system, for backrest section lowering;
- Transfusion / infusion armrests;
- Lateral support with instrument tray holder;
- Colposcope arm support, with anti-crushing safety stop;
- No. 2 additional electric sockets in the base (230 V - 300 VA);
- Foot control cable wind up;
- Viscoelastic polyurethane upholsteries (memory-foam), free of seams;
- Liquid containment basin with valve, for urological check-up;
- Electrically warmed seat upholstery;
- Back-up battery (for brief movements in absence of power supply).

The operation and the performance of these accessories are described in chapter 7.

2.3. WARRANTY.

NOTE.

Upholstery damage caused by contact with clothes or clothing accessories (ex: belts) having metal, sharp, pointed or similar parts will not be replaced under warranty.

CAUTION.

The use of disinfectants or detergents with a strong alcoholic content on the plastic components of the equipment may damage the surface of the product, especially the upholstery.

2.4. IDENTIFICATION AND LABELLING.

The identification label is positioned externally on the chair. It indicates the chair model, serial number and other information. Figure 1 highlights the label position and a facsimile.


FIG. 1

	Manufacturer		Manufacturing date
	Model. It identifies unambiguously the device of the product family		Serial number of the device
	Medical device		UDI code. UDI coding may not be shown on the label
	Read the operating instruction carefully		Contains type B applied parts
	CE marking		Dispose of the product in compliance with the procedures provided for electric and electronic equipment

2.5. TECHNICAL SPECIFICATIONS.


2.5.1. Sector-specific standard compliance.

The gynecological patient chair model SERENITY NEXT complies with with the Regulation (EU) 2017/745 (MDR), concerning medical devices.

In particular, design and manufacture of the equipment comply with the following international standards:

- EN 60601-1:2006 / A11 A1;
- EN 60601-1-2:2015.

2.5.2. Technical data.

Power supply – voltage	V	230 AC
Power supply – frequency	Hz	50 – 60
Input power (including auxiliary sockets)	VA	510 (1100)
Maximum output power available for auxiliary socket	VA	300
Main fuses		2 x T 3,15A; 250V
Controls voltage	V	5
External controls interface		Yes
Insulation class		I
Applied parts type		B
Medical device risk class		I
Maximum height	mm	920 (1095)
Minimum height	mm	600 (545)
Width (without accessories)	mm	750
Length (with leg rest completely extracted)	mm	1840
Backrest inclination	°	0° – 65°
Trendelenburg inclination (backrest and seat)	°	-10° - +15°
Maximum basin capacity	l	2,5
Gross weight	Kg	180 (depending on the accessories included)
Net weight	Kg	150 (depending on the accessories included)
Packaging dimensions	mm	1430x810x990 (h)
Maximum lifting capability (patient + unit)	Kg	250
Operation mode		Intermittent loading, 2 min. ON – 18 min. OFF
Minimum space required for installation	m	3 x 2

2.5.3. Electromagnetic compatibility.

The SERENITY NEXT chair does not cause influences of an electromagnetic nature or other type with other equipment, nor is it susceptible to influences from other devices, and is compliant with the requirements of the harmonized standard EN 60601-1-2. The following tables provide details on the electromagnetic environment in which the product can be used.

Emission aspects.

Emission test	Conformity	Electromagnetic environment - guide
RF Emissions Cispr 11	Group 1	The product uses RF energy for its internal function only. Therefore, its RF emissions are extremely low and not likely to cause interference with nearby electronic devices.
RF Emissions Cispr 11	Class B	It is possible to use the equipment in all buildings, including domestic use buildings, and those directly connected to a public power supply network with a low voltage supplying buildings for domestic use, taking installation steps, and guaranteeing an increased installation distance for potentially sensitive devices.
Harmonic emissions IEC 61000-3-2	Class A Conformity	It is possible to use the equipment in all buildings, including domestic buildings, and those directly connected to a public power supply network with a low voltage supplying buildings for domestic use.
Emission fluctuations voltage/flicker IEC61000-3-3	Conformity	

Immunity aspects.

The SERENITY NEXT chair is intended for use in an electromagnetic environment under specification. The customer or user should ensure that it is used correctly in such an environment.

Immunity test	Test level EN 60601-1-2	Conformity level	Electromagnetic environment - guide
Electrostatic discharge (ESD) EN 61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If the floors are covered with synthetic material, the relative humidity should be at least 30%.
Burst / Fast transient EN 61000-4-4	±2kV power supply lines	±2kV power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode ±2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	< 5% UT (>95% dip UT) for 0,5 cycles 40% UT (60% dip UT) for 5 cycles 70% UT (30% dip UT) for 25 cycles < 5% UT (>95% dip UT) for 5 seconds	< 5% UT (>95% dip UT) for 0,5 cycles 40% UT (60% dip UT) for 5 cycles 70% UT (30% dip UT) for 25 cycles < 5% UT (>95% dip UT) for 5 seconds	Mains power quality should be that of typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency magnetic field EN 61000-4-8	3 A/m	3 A/m	Magnetic power frequency fields should be that of a typical commercial or hospital environment.

RF immunity aspects.

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an electromagnetic environment.

Immunity test	Test level EN 60601-1-2	Conformity level	Electromagnetic environment - guide
RF Conducted EN 61000-4-6	3 Veff from 150kHz to 80MHz	3 Veff from 150kHz to 80MHz	Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distances calculated from the equation, applicable to the frequency of the transmitter. Recommended separation distance: $d = 1,2 \cdot \sqrt{P}$ from 150kHz to 80MHz $d = 1,2 \cdot \sqrt{P}$ from 80 MHz to 800 MHz $d = 2,3 \cdot \sqrt{P}$ from 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).
RF Radiated EN 61000-4-3	3 Veff from 80MHz to 2,5GHz	3 Veff from 80MHz to 2,5GHz	
Field strength from fixed RF transmitters, as determined in an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:			

Recommended separation distance between portable and mobile RF communication equipment and the device.

The SERENITY NEXT chair is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the device as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of the transmitter (m)		
	From 150kHz to 80MHz $d = 1,2 \cdot \sqrt{P}$	From 80MHz to 800MHz $d = 1,2 \cdot \sqrt{P}$	From 800MHz to 2GHz $d = 2,3 \cdot \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters with a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum power rating of the transmitter in Watts (W) according to the transmitter manufacturer. Note: (1) At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. (2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

3. INSTALLATION

3.1. UNPACKING.

To remove the cardboard box, it is necessary to cut the straps and to unscrew the screws located in the lower part of the box.


FIG. 3

3.2. SUPPLY STATUS.

The following elements must be contained inside the packaging:

- A SERENITY NEXT chair;
- A paper towel roll;
- An operating instruction manual;
- A bottle of fabric cleaner.

In case of mistake or missing parts, contact the dealer to ask for assistance.

3.3. HANDLING.

To remove the chair from the wooden pallet, it is necessary to remove the two screws which fix the chair base to the pallet, situated in the front side of the base (see figure 3).

To move the chair from the wooden pallet, one or more assistants are needed.

In its front part, the chair can be displaced by grabbing the side instrument holders or the operator's foot rest (G).

For the rear part, a transport handle (F) is available. It can be assembled between the two lifting columns and has to be re-moved once the chair is displaced.

3.4. INSTALLATION.

3.4.1. Typical installation.

Proceed as follows:

1. Be sure that the installation site conforms to local regulatory technical requirements for medical studios or outpatient clinics;
2. Be sure that the electrical system of the cabinet has a safe earthing and a residual-current device system;
3. Be sure that the main supply of the electrical wiring of the cabinet corresponds to the one of the identification label;
4. Check that the fuses, located in the lower part of the base, are correctly fixed, because during transport they could have become loose;
5. Connect the power cable firstly to the chair and then to the power socket;
6. Activate the main switch to turn on the chair. Now the chair is ready for use. In case of non operation of the chair, refer to section 6 about troubleshooting;
7. If needed, drive the chair to a comfortable height and assemble, if provided, the accessories.

NOTE.

Equipment must not be used in potentially explosive environment, or in rooms not protected against weather conditions. Not suitable for use in altered atmospheres.

3.4.2. Position choice and stability.

To ensure a good chair installation, please verify the stability and flatness of the floor in the location where the product will be installed: if necessary, a damaged or unstable floor should be repaired.

The chair is designed to remain stable without floor fixing.

To reach that aim, it is provided with four adjustable feet (figure 4) to level the base in case of slightly odd surfaces.

The product is designed for stationary installation.

Upon request at the moment of the order, as an alternative to the levelling feet is however possible to equip the chair with braked castors, which allows short range movements. For further detail refer to chapter 7.

NOTE.

The levelling operation must be conducted during the installation phase.


FIG. 4

WARNING.

Tecnodent declines all and any responsibility for damage caused by the non-compliance with the above instructions.

4. USING THE PRODUCT

The advanced technology applied to the production of this apparatus has allowed the design of a patient chair with multiple functions and performances, all of them easily usable: all regulations and functions of the chair are coordinated, processed and stored in a "data bank", by a microprocessor inside the chair.

WARNING.

Before using the information contained in this section, make sure that the chair is correctly installed, by the service technician.

CAUTION.

Use in a clean place (specialized surgery, clinic) and avoid accidental spilling or falling of liquid substances on the product.

CAUTION.

In case of a product malfunction that is not described in this documentation, put equipment out of service (shut it down, unplug main supply, mark the equipment as "out of order" or similar) and call the technical service.

NOTE.


The chair is provided with an electronical system which detects a possible overload of the motors. This system is designed to intervene in situations like excessive weight on the chair or mechanical locks of the movements.

4.1. TYPICAL USE.

To use the product, the operator must be in the proximity of the chair, while acting on the controls. Care should be taken throughout the operation of the product, as the presence of automatic movements could cause hazardous situations, even if moderate level ones.

The chair is provided with several emergency stop systems, e.g., in case of presence of obstacles during the movements, or if the operator checks any unsuitable situation.

4.2. CONTROLS.

The controls are shown in figure 5.

- P Main switch (Power on/off).
- L Push-button hand-held panel. It controls the raising, lowering and backrest tilting, and Trendelenburg movements.

4.3. MOVEMENTS.

- To raise the backrest, keep pressed the left part of the push-button A.
- To lean the backrest, keep pressed the right part of the push-button A.
- To raise the motorised leg supports (if present), keep pressed the left part of the push-button B.
- To lower the motorised leg supports (if present), keep pressed the right part of the push-button B.
- To raise the chair, keep pressed the left part of the push-button C.
- To lower the chair, keep pressed the right part of the push-button C.
- To lean the chair towards the Trendelenburg position, keep pressed the left part of the push-button D.
- To make the chair return from the Trendelenburg position towards the patient loading assistance position (anti-Trendelenburg), keep pressed the right part of the push-button D.

5. MAINTENANCE

5.1. GENERAL INFORMATION.

SERENITY NEXT chairs were conceived to not require any servicing during their working life. Consequently, it is not required to perform any adjustment of the equipment. If the chair will not be used for several days, it is recommended to cover it with a cloth, to turn off the main switch, and to disconnect the plug from the socket.

If the SERENITY NEXT is used according to the technical specifications and operating instructions specified in this manual, the expected service life has been calculated in at least 10 years, based on laboratory simulations.

After that time, the chair needs to be serviced by qualified personnel, especially the electronic parts and the ones subject to wear.

WARNING.

Qualified service technicians must carry out any technical adjustments or repairs not listed in this chapter.

5.2. CLEANING.

For a better and longer life of the equipment, it is necessary to perform an accurate, methodical and periodical general cleaning of the chair. It is recommended to proceed as follows.

CAUTION.

The use of disinfectants or detergents with a strong alcoholic content on the plastic components of the apparatus may damage the surface of the product, especially the upholstery.
It is recommended to avoid the use of any dissolvent, strong detergent or abrasive agent, even to remove "difficult" stains.

5.2.1. Upholstery.

The fabric that covers the chair upholstery has to be cleaned using the liquid contained in the bottle kit supplied with the chair, using a soft cloth to avoid any surface scratching and to guarantee a better elasticity and smoothness of the surface.

NOTE.

Once the liquid detergent has been consumed, please contact the nearest assistance centre for new product.


5.2.2. Plastic or metallic parts.

The polyurethane parts, as well as the metal parts (painted and not) have to be cleaned with a soft cloth dipped in water and soap.

5.3. MAIN FUSES REPLACEMENT.

The SERENITY NEXT chair is equipped with a double overloading protection.

The two fuses are located in the lower part of the base, over the main switch, as indicated in figure 6.

To replace the fuses, proceed as follows:

1. Turn off the chair's main switch;
2. Disconnect the plug from the mains socket;
3. Unscrew the fuse holder cap, by using a crosshead screwdriver of medium size;
4. Replace the damaged fuses;
5. Reassemble the components, following the inverse sequence of the above steps.

WARNING.

Fuses must only be replaced by other fuses having the same nominal characteristics. Table "Technical data" in paragraph 2.5.2 contains the relative specifications.

5.4. UPHOLSTERY REPLACEMENT.

To replace the chair upholsteries, proceed as follows and as shown in figure 7.


FIG. 7

5.4.1. Seat.

The seat upholstery is press fitted by means of four plastic press-fit plugs. To remove simply raise the upholstered section applying a light force until it detaches from the corresponding mated component.

To attach, repeat the following procedure in reverse order.

WARNING.

If the chair is provided with the electrically warmed seat upholstery, please refer to paragraph 7.6.

5.4.2. Backrest.

The seat upholstery is press fitted by means of four plastic press-fit plugs. To remove simply raise the upholstered section applying a light force until it detaches from the corresponding mated component.

In order to re-assemble it, lay the upholstery to the back rest frame, align the snap-on caps with the relevant pins, and press in order to hook them in.

6. TROUBLESHOOTING

The steps to detect a faulty condition on the SERENITY NEXT chair are described in table-form. Each table reports the title of the troubleshooting procedure, and three columns whose content can help to detect and solve the problem. Below is an example:

Procedure Name		
Problem	Test Procedure	Corrective Action(s)
Patient chair problem description.	If the answer to the proposed question is "YES" look at the question written in the next box below; if the answer to the question is "NO" take all the corrective action(s) described in the right box before moving to the next question.	Action(s) to be carried out in order to fulfil the requirements of the question asked by the test procedure at the side. Only after all actions have been carried out with positive result it is possible to move to the next test procedure.

The table-form takes into consideration only electrical or electronic faulty conditions as mechanical faults are easily detectable by a simple visual inspection. In order to execute a visual inspection, it is sufficient to remove the upholstery parts, as described in paragraph 5.4.

Main supply control procedure		
No chair movement.	Is the supply voltage correctly rated?	If an OHM-meter or a similar testing device is not available, plug any other apparatus (i.e. a lamp) into the socket and verify the presence of the power supply.
	Is the power cord plugged in?	Plug the cord into the socket.
	Is the green light of the main switch on?	Turn the main switch on.
		After having unplugged the chair, check the main fuses. If necessary, replace them as described in the user manual. Otherwise call the technical service.
	Does the PCB emit a double "BEEP" sound when the main switch is turned on?	The PCB has probably failed. Call the technical service.
	The problem is not related to the chair main supply group. Continue checking the chair movements one by one.	

Chair lifting movement control procedure		
The chair does not lift.	Is there too much weight on the chair?	Apply the chair movement testing procedure.
	Remove the excessive weight. Look at the chair operating instruction manual to determine the maximum chair lifting capability.	

Backrest tilting movement control procedure		
The chair backrest does not tilt.	Is there an object that collides with the chair, activating the safety device and stopping the backrest tilt movement?	Apply the chair movement testing procedure.
	Remove the object that collides with the chair.	

Backrest raising movement control procedure		
The chair backrest does not raise.	Was the chair movement testing procedure duly applied?	Apply the chair movement testing procedure.
	Call the technical service.	

Chair movement testing procedure		
The chair does not move at all, even if the main switch is switched on and the mains group is working correctly.	Do you hear a "CLICK" sound when you try to activate the controls?	If no "CLICK" sound is emitted probably the PCB has failed. Call the servicing technicians.
	Does the motor hum?	The motor may have over-heated and the thermal security switch may have turned on. Wait for about 15 minutes, allowing the motor to cool down and operate the motor again. If the chair does not move or the motor overheats immediately after, call the servicing technicians.
	It is impossible to determine the cause of the failure unless further and more complex testing is applied. Call the servicing technicians.	

Chair programming control procedure		
The chair moves properly in all manual movements, but it is impossible to store working position programs.	Was the programming procedure properly applied, as explained in paragraph 7.1.2?	Apply the procedure as quoted in paragraph 7.1.2.
	Did you hear a "BEEP" sound when the programming procedure was completed, to confirm the programs' saving?	It is impossible to determine the cause of the failure unless further and more complex tests are performed. Call the technical service.
	The PCB has probably failed. Call the technical service.	

NOTE.


If the foot control display shows an error message ("ERR" followed by a number code), contact the technical assistance, referencing the visualised code.

7. ACCESSORIES

7.1. FOOT CONTROL WITH DISPLAY.

The SERENITY NEXT chairs can be equipped, upon request while ordering, with a foot control (shown in figure 8) that allows a more advanced use of the product. In addition to let the professional have his hands free and to show possible error messages, it allows multiple functions:

- Automatic return to the patient exit position ("zero");
- Programmable working positions;
- Work profiles;
- "SLOW MODE" setting.

7.1.1. Basic movements.

RAISING AND LOWERING.

The symbol on the left identifies the movement.
Keep the lower key pressed to lower the chair.
Keep the upper key pressed to raise the chair.

BACKREST INCLINATION.

The symbol on the left identifies the movement.
Keep the lower key pressed to raise the back section.
Keep the upper key pressed to lower the back section.

TRENDELENBURG AND LOADING ASSISTANCE (ANTI-TRENDELENBURG).

The symbol on the left identifies the movement.
Keep the lower key pressed to incline the upper part of the chair forwards.
Keep the upper key pressed to incline the upper part of the chair backwards.

PROGRAMS.

Tap one of the keys P1 or P2 or P3 and the chair will automatically return to the programmed working position.
Tap the Ø key and the chair will automatically return to the patient exit position.

NOTE.


To distinguish between "keep pressed" and "tap", the controller needs a delay of approximately 0.5 seconds.

NOTE.

Any switch, if pushed during an automatic movement, works like a safety stop.

7.1.2. Programming.

1. Set the chair to the patient exit position by tapping the Ø button on the foot control.
2. Configure the chair to the desired position by activating the corresponding foot controls (see paragraph 7.1.1).
3. Press the key shown on the left and, within 8 seconds, press the key "P1" or "P2" or "P3".

An acoustic signal ("BEEP") indicates that the work position has been memorized.

Every time that the controls "P1", "P2" or "P3" are tapped, the chair will automatically configure to the pre-programmed working position.

NOTE.

To avoid losing the programmed position, it is recommended to regularly set the chair to its "zero" position, by using the "Ø" function at any patient change.

7.1.3. Work profiles.

The key indicated in the figure enables the selection of three work profiles.

Tap the key to switch between three work profiles: for each profile it is possible to memorize three work positions (corresponding to the keys P1, P2 and P3). Pressing the key permits access to the successive profiles before eventually returning to the first one.

The foot control display shows the selected profiles (USER1, USER2, USER3), while the transition of each profile is indicated by a single acoustic signal ("BEEP").

7.1.4. "SLOW MODE".

Keeping the key pressed for approximately two seconds, enters the chair into SLOW MODE, as indicated by the foot control display. While this setting is active, all the movements are significantly slower; e.g. to assist specific therapies.

To exit from this setting, keep again the key pressed. The activation and deactivation are indicated by four quick intermittent acoustic signals ("BEEP-BEEP-BEEP-BEEP").

7.2. PROTECTIVE PAPER TOWEL HOLDER.

The chair is fitted with a system to house rolls of protective paper, essential for patient hygiene, forming a partial layer of protection between the patient and device.

The protective paper is placed in the upper part of the back section and covers the entire chair upholstery.

NOTE.

For improved hygiene between patients, dispose of the protective paper used during the examination and withdraw fresh sheets from the roll.

Once the paper roll is finished, it can be replaced without specific procedures: remove the pin, replace the paper towel roll and place again the pin in its housing.

7.3. UPHOLSTERY WITH FIRE RETARDANT AND ANTI-STATIC LEATHERETTE.

Upon request at the moment of the order, the chair can be produced with a specific fire-retardant leatherette which can allow to obtain the fireproof class 1/IM according to the standard UNI 9175. This option can be retrofitted after the purchase, by ordering the specific spare part kit.

7.4. MOULDED MEMORY FOAM UPHOLSTERY

Upon request while ordering, the upholstery can be produced in viscoelastic polyurethane memory foam, without seams.

This option can be retrofitted after the purchase, by ordering the specific spare part kit.

This upholstery offers greater comfort and adapts better to the patient's body thanks to the "memory" effect, the slow return of its elasticity and the fact that its deformability depends on body heat.

Moreover, the anti-decubitus properties of the material make it ideal for lengthy treatments.

Finally, the moulding process and the lack of stitching, a unique feature of Tecnodent memory foam, makes it possible to maintain an excellent level of hygiene.

7.5. LIQUID CONTAINMENT BASIN, WITH VALVE.

If specified while ordering, the chair can be equipped with a basin which allows the containments of liquids, with a valve. This containment basin can be connected to a drain system. This accessory allows the use of the chair for urological examinations.


FIG. 9

7.6. ELECTRICALLY WARMED SEAT UPHOLSTERY.

Upon request while ordering, the chair can be equipped with a heating system in the seat upholstery, allowing the set-up with three temperature levels.

OPERATION.

This accessory can be controlled by a dedicated hand-held control. On this control, four buttons are available; by pressing those buttons it is possible to adjust the power of the seat warming.

MAINTENANCE.

The presence of this accessory needs some additional operations if it is necessary to remove the seat upholstery.

- Turn completely off the chair, including the back-up battery, if provided.
- Raise the seat upholstery, applying a light force until it detaches from the corresponding mated component, but do not remove it yet.
- Unplug the power supply connector.
- Remove the upholstery.

To reassemble the part, follow the above instructions in reverse sequence.


FIG. 10

7.7. BACK-UP BATTERY.

Upon request while ordering, the chair can be equipped with an internal back-up battery.

This accessory allows to use the basic movements of the chair, even in absence of power supply. For example, the back-up battery, in case of power supply failure, allows to drive the chair into a position from which the patient can be unloaded safely.

NOTE.


In absence of power supply, the chair disables all the advanced functions, such as user-defined programs, user profiles, Trendelenburg movement, warmed seat upholstery.

OPERATION.

Should the power supply be interrupted, including by simply turning the mains switch off, the foot control display shows the message BATTERY.

CHAIR POWER OFF.

In order to completely turn off the chair, it is needed to keep pressed the MEMORY button for five seconds. Without this procedure, the chair will completely turn off autonomously after five minutes of chair inactivity (no movement operated).

NOTE.

Once the chair has been completely turned off, it is only possible to turn it on when the power supply is restored, by activating the main switch.

7.8. MAGNETIC CUSHION.

The SERENITY NEXT chairs may be fitted upon request with a magnetic cushion (see figure 11) for improved patient comfort.

MAINTENANCE.

The head support upholstery may be easily replaced by separating it from its counterpart.



FIG. 11

7.9. INTEGRATED RETRACTABLE LEG REST.

This accessory is designed in order to allow, during a visit, procedures which require the elongation of the patient legs. It allows to temporarily transform the chair into an examination bed.

The presence of this accessory can be specified only at the moment of the order.

If present, the leg rest section is housed under the seat.

UNFOLDING.

- Grab the front border of the leg rest and extract it by completely extending its guides (1).
- Once the guides are completely extended, raise slightly the front border and continue the movement (2) in order to toggle the lifting mechanism.
- Press the lever on the leg rest guides (L) and at the same time push gently the leg rest to the seat section, until the first mechanical stop is reached, near its final position (3).

FOLDING.

Once the examination is concluded, and if there is any necessity to unload the patient, it is necessary to place the leg rest back into its housing:

- Move gently the leg rest away from the seat, until a locking sound is heard (4).
- Drive gently the leg rest (5) to its housing under the seat, and fold the mechanism until a locking sound is heard.

CAUTION.

The following operations could damage the equipment:

- Overload the leg rest (e.g., by a person sitting on the front side);
- Let the patient get on or off while the leg rest is horizontal;
- Activate the chair lowering movement while an obstacle is under the leg rest section.

WARNING.

Keep clear the area in front of the chair, under the leg rest.

7.10. PATIENT SIDE HANDLES

Upon request, while ordering, the chair may be fitted with a pair of handles, that have a double function:

- To allow the patient to hold better during the chair movements;
- To house the transfusion armrest (see paragraph 7.12).

The handles could be rotated downwards, in order to allow an easier access to patient area if the operator is near the side of the chair.

To lower the handles, pull them outward and rotate them by 180°.

NOTE.

Patient handles are not available if the chair is provided with motorised leg rest supports.

7.11. SIDE INSTRUMENTS HOLDER.

Upon request, at the moment of the order, the chair may be fitted with a pair of lateral supports, designed to house various instruments.

7.12. TRANSFUSION / INFUSION ARMREST.

The chair may be fitted with an armrest specifically designed for blood transfusion or infusion. The component must be mounted onto one of the two handles. To attach it, simply hook it to the handle and tighten the knob A as shown in fig 13.


FIG. 13

7.13. LEG SUPPORTS.

To hold and move the patient's legs, two kinds of supports can be provided: a model is equipped with under-thigh holders and a model is equipped with foot holders. The housings for these supports are on the side of the chair, under the seat upholstery.

7.13.1. Under thigh (Göpel) model.

Figure 14 shows the support with under-thigh holder. To adjust the inclination of the holder (1), loosen the handle (2), lead the holder to the desired position, then tighten again the handle (2). This model of leg holder could be provided with different kinds of support, depending on the operator's needs.

TELESCOPIC SUPPORT (WITH 3 ADJUSTMENTS).

To assemble these supports on the chair, it is necessary to loosen the handle on the housing A in the chair, insert the knurled tail in this housing, drive the support to the desired position, then tighten the handle again.

In order to adjust the supports, proceed as follows:

1. Raise the patient leg;
2. Loosen the handle (3) and modify the support length;
3. Once the desired height is attained, tighten the lever (3) again.


FIG. 14

VERTICAL FOLDABLE SUPPORTS (WITH 2 ADJUSTMENTS).

This accessory is advised in case the operator needs to access the patient frequently from a sideways position and does not intend (or has no space) to remove the telescopic supports.

To adjust the height, loosen the handle on the chair housing (A), reach the desired position and tighten the lever again.


FIG. 15

7.13.2. Under feet model.

Figure 16 shows the leg holder with support under the foot. This accessory can be provided only with telescopic support. To adjust the holder, proceed as follows:

1. Loosen the handle (1) and modify the support length;
2. Once the desired height is attained, tighten the handle (1) again;
3. To adjust the holder inclination (2), loosen the handle (3), drive the part to the desired position, then tighten the handle (3) again.


FIG. 16

7.13.3. Motorised leg rest supports.

Upon request at the moment of the order, the chair can be equipped with leg holders with under feet support and motorised movement. Two kind of movements are available: lifting / lowering (electrically assisted) and rotation (manual).

OPERATION.

The leg rests can be controlled by a keypad with two buttons, placed on the side of the seat, and with a handle placed on the rotation pin.

- To lift the leg holders, press the corresponding button on the keypad.
- To lower the leg holders, press the corresponding button on the keypad.
- To rotate the leg holders, loosen the lever, manually rotate the holder to the desired position, then tighten the lever again.


FIG. 17

7.14. INSTRUMENT TRAY HOLDER SUPPORT.

Figure 18 shows the instrument tray holder.

This accessory is fitted to the colposcope arm tube, fixed with the provided screws A and B. Depending on the accessories assembled on the chair, specific adaptors may be required; they can be provided upon request.


FIG. 18

7.15. COLPOSCOPE ARM SUPPORT.

To make possible the attachment of a colposcope to the chair, Tecnodent can provide a specific support, shown in figure 19. The accessory must be inserted into the housings placed under the seat frame, on the front of the chair, then fixed in place by tightening one of the ends with the provided screw. Depending on the accessories assembled on the chair, specific adaptors may be required; they can be provided upon request.


FIG. 19

7.16. “QUICK RELEASE” SYSTEM FOR BACKREST LOWERING.

If specified while ordering, the SERENITY NEXT chair may be fitted with a mechanism that allows a quick lowering of the backrest section. The control knob is placed on both sides of the backrest. To unlock the system, rotate the knob as shown in figure 20, while holding the weight of the backrest.

Once the need of quick backrest lowering has ended, it is necessary to perform a complete backrest movement by using the electrical controls of the chair, aimed to bring the backrest itself to its end positions. Such precaution will prevent the possibility of failures, especially with pre-programmed working positions.


FIG. 20

7.17. ADDITIONAL ELECTRIC SOCKET ON THE BASE.

Upon request, at the moment of the order, Tecnodent can equip the device with two additional sockets, located at the back of the chair base, besides the mains switch.

The protection fuses are placed near the socket themselves.

To replace the fuses, operate in a similar way, as specified for the chair power supply fuses (paragraph 5.3).

WARNING.

The fuses protect the overloading of the additional sockets. It is recommended to verify that the power absorption of the equipment connected does not exceed the maximum power rated for the socket (300 VA).


FIG. 21

7.18. FOOT CONTROL CABLE WIND UP.

Upon request, at the moment of the order, the chair may be fitted with a foot control cable wind up mechanism. This device allows enhanced order in the zone where the foot control is placed, by recalling the cable inside the base.


FIG. 22

7.19. BASE WITH WHEELS.

As an alternative to the adjustable feet, upon request while ordering, the chair may be fitted with braked castors, as shown in figure 23.

CAUTION.

This accessory allows only short-range movements, for example during cleaning operations.



FIG. 23

The brake of the castors can be activated or released by using the levers P, available on each wheel.

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Tecnodent S.r.l.