



Zerlegbare Rongeurs mit Silikongriff

Dismantable Rongeurs with Silicone Handle

Pinzas saca-bocadas desmontables con mango de silicona

Rongeurs démontables avec manche en silicone


Pinze per laminectomia smontabili con manico in silicone



Medicon ...DieBLAUEN.2



Zerlegbarkeit gewährleistet optimalen Zugang für die Aufbereitung
Dismantling ensures optimal access for reprocessing
Desmontaje que garantiza un acceso óptimo en el reprocesamiento
Démontage garantissant un accès optimal pour le traitement
La facilità di smontaggio assicura l'accesso ottimale per il ricondizionamento




- ▶ Einteilige Zerlegetechnik
- ▶ Ergonomisches Design
- ▶ Silikongriff für höchsten Arbeitskomfort
- ▶ One piece dismantable rongeur
- ▶ Ergonomic design
- ▶ Silicone handle for high working comfort
- ▶ Pinza saca-bocada desmontable de una sola pieza
- ▶ Diseño ergonómico
- ▶ Mango de silicona para la máxima confort de trabajo
- ▶ Technique de démontage en une pièce
- ▶ Design ergonomique
- ▶ Manche en silicone pour un confort de travail maximal
- ▶ Tecnica di smontaggio monopezzo
- ▶ Design ergonomici
- ▶ Manico in silicone per altissima comfort lavorativo



Lange Standzeit, hohe Oberflächenhärte und blendfreie Oberfläche
Long lasting cutting ability, high surface hardness and non-dazzling surface
Máxima duración del corte, alta dureza de la superficie y superficie antideslumbrante
Longue durée de vie, grande dureté de la surface et surface non-éblouissante
Massima durata del taglio, alta durezza della superficie e superficie antiriflesso

Zerlegbare Rongeurs mit Silikongriff
Dismantable Rongeurs with Silicone Handle
Pinzas saca-bocadas desmontables con mango de silicona
Rongeurs démontables avec manche en silicone
Pinze per laminectomia smontabili con manico in silicone



Schaftlänge Working length Longitud de trabajo Lunghezza di lavoro	Maulbreite Jaw width Ancho de la boca Largeur des mors Larghezza del morso	Maullänge Jaw length Longitud de la boca Longueur des mors Lunghezza del morso	gerade straight recta droit rett	aufgewinkelt angled upwards angulada hacia arriba coudé en haut angolata verso l'alto	abgewinkelt angled downwards angulada hacia abajo coudé en bas angolata verso il basso
a	b	c	0°	30°	30°

	130 mm	2 mm	10 mm	57.37.10	57.37.12	57.37.14	Cushing
	155 mm	2 mm	5 mm	57.37.00	57.37.01		Micro Cushing
	155 mm	2 mm	5 mm	57.37.02	57.37.05		Micro Cushing
	155 mm	3 mm	10 mm	57.37.30			Love Gruenwald
	155 mm	2 mm	12 mm	57.37.55			Caspar
	155 mm	3 mm	12 mm	57.37.56			Caspar
	155 mm	2 mm	14 mm	57.37.50	57.37.53		Caspar
	155 mm	3 mm	14 mm	57.37.51	57.37.54		Caspar
	155 mm	4 mm	14 mm	57.37.52			Caspar
	155 mm	4 mm	14 mm	57.37.57			Caspar
	180 mm	2 mm	10 mm	57.37.03	57.37.04		Cushing
	180 mm	2 mm	10 mm	57.37.11	57.37.13	57.37.15	Cushing
	180 mm	3 mm	10 mm	57.37.31	57.37.32	57.37.33	Love Gruenwald
	180 mm	4 mm	10 mm	57.37.40	57.37.41	57.37.42	Spurling
	180 mm	2 mm	12 mm	57.37.58			Caspar
	180 mm	3 mm	12 mm	57.37.59	57.37.61		Caspar
	180 mm	4 mm	14 mm	57.37.60	57.37.62		Caspar

magic cut

magic cut



Maulteile gefenestert und gezahnt
 Jaws fenestrated and serrated
 Bocas fenestradas y dentadas
 Mors fenêtrés et dentés
 Morsi fenestrati e dentellati



Maulteile gezahnt
 Jaws serrated
 Bocas dentadas
 Mors dentés
 Morsi dentellati



||| CMF

||| SURGICAL

||| NEURO+SPINE

Medicon eG

Gänsäcker 15
D-78532 Tuttlingen

P. O. Box 44 55
D-78509 Tuttlingen

Tel.: +49 (0) 74 62 / 20 09-0
Fax: +49 (0) 74 62 / 20 09-50

E-Mail: sales@medicon.de
Internet: www.medicon.de

Germany

Not all products illustrated in this documentation are available for sale in all countries. Please contact your local distributor for more information concerning availability and regulatory information (e.g. classification) of the mentioned products.



Gedruckt in Deutschland
Printed in Germany
Impreso en Alemania
Imprimé en Allemagne
Stampato in R.F.G.
Trykkes i Tyskland

Modelländerungen vorbehalten
Patterns are subject to change
Salvo modificaciones
Tous droits réservés des changements de modèle
Ci riserviamo la facoltà di cambiamenti nei modelli

451.04.58

© Copyright 06/2024, MEDICON eG, Tuttlingen, Germany

medicon[®]

||| NEURO+SPINE

