HEINE headlights

Overview

It goes without saying that all of our headlights are premium products. But, still, they differ in one or two details. Read on for a summary of what these different details are and all their benefits:





	MicroLight2	ML4 LED HeadLight
Illumination	LEDHO	
Light intensity*	typ. 55,000 lux	typ. 90,000 lux
Colour temperature (kelvin)	typ. 4,000 K	typ. 5,500 K
Spot size	81 mm	30 – 80 mm
Power sources operating time	mPack mini typ. 9 hours	mPack UNPLUGGED typ. 3.5 hours mPack typ. 8.5 Std.
Weight	80 g with S-FRAME spectacle frame 190 g with Lightweight headband	260 g 460 g with mPack UNPLUGGED
Angle of declination	fully-adjustable	
Brightness control	0 – 100 %	
Features	waterproof, optional yellow filter, head- or spectacle-worn,	adjustable spot size, optional yellow filter, polarisation filter, expandable with HR and HRP loupes, head-worn,
Catalogue page	139-140	143-144

 $^{^*\}mbox{Measured}$ at a working distance of 180 mm from the light emission to the object.




HEINE ML4 LED HeadLight Ultra bright coaxial LED^{HQ} illumination



- ► There's LED. And there's LED^{HQ}. HEINE has set a new standard which maintains that only the best is good enough, from the selection of materials to processing, from light intensity to dimmability, and from thermal management to a colour rendering index as high as possible. That's what we call LED in HEINE quality or LED^{HQ}
- ► Adjustable illumination spot size: 30 mm to 80 mm spot size range (at 420 mm working distance) to fit any examination situation
- Bright and homogenous illumination. Absolutely bright light spot that is uniform from edge-to-edge for the perfect illumination in all examination situations
- ► Flexibility. Flexible power source options: 100 % mobility with a choice of cable-free mPack UNPLUGGED head-worn battery or mPack belt-worn battery pack
- ► Comfortable and secure fit. The Professional L headband has multiple adjustment points and soft padding for comfort and stability even during long examinations
- ▶ **Optional filters.** Polarisation filter P2 improves contrast and yellow filter reduces blue light

Technical specification



[01] Competitor product [02] HEINE ML4 LED HeadLight

typ. 90,000 lux*
typ. 5,500 kelvin
typ. 50,000 hours
30mm to 80mm at 420mm working distance
stepless and individually adjustable for every examination
headband-mounted
typ. 8.5 hours with mPack and typ. 3.5 hours with mPack UNPLUGGED of continuous ON time at full power

^{*}Measured at a working distance of 180 mm from the light emission to the object.

HEINE ML4 LED HeadLight

Power options

HEINE mPACK



- ▶ Portable power pack. Does not tie the user to desk, wall or static power
- ▶ Li-ion technology. Triple capacity compared with current NiCd systems. No memory effect, autonomy for typ. 8.5 hours of continuous ON time at full power
- Fast charge. Charges to full capacity in typ. 2 hours
- ► Charge status indicator. Always lets you know how much operating charge you have left
- Brightness control. Optimal adjustment of light intensity
- Flexible charging options. EN 50 desk and wall charger, or plug-in trans-
- ▶ Intelligent. Can be used as mains power supply even while in "charge" mode
- 100 240 V compatible transformer incl. 4 plugs for world wide compatibility

HEINE mPACK UNPLUGGED HEADBAND-MOUNTED RECHARGEABLE BATTERY



- ▶ Wireless. 100 % mobility for absolute freedom of movement
- Rechargeable battery integrated on the headband
- Lightweight with balanced weight distribution. High level of comfort
- Charge status indicator indicates remaining operating time
- Intelligent. Can be used as mains power supply even in "charge" mode
- Fast charge. Charges to full capacity in typ. 2 hours. Operating time typ. 3.5 hours of continuous ON time at full power
- Top safety standards. Integrated charge and performance checks during operation. Lithium polymer technology. UL compliant










HEINE ML4 LED HeadLight

Ordering information



$\ensuremath{\mathsf{ML4}}$ LED HeadLight with Professional L headband





ML4 LED HeadLight accessories	
detachable yellow filter [06]	J-000.31.321
extension cord from plug-in transformer UNPLUGGED to mPack UNPLUGGED, 2 m [07]	X-000.99.668
polarisation filter P2 for ML4 LED and HR loupes [08] + [09]	J-000.31.326
spare polarisation filter P2 for HR loupes (2 pairs) [09]	C-000.32.536

[03]

Power supplies may deviate in size and shape.











